

DNA CREATING A CODIS CMF FILE

- A.1. Open the Genemapper *ID-X* project.
- A.2. Select the **SAMPLES LINK** in the **ANALYSIS SUMMARY** tab.
- A.3. Select the **CODIS EXPORT TABLE SETTING** from the drop down list.
- A.4. For each row select the appropriate setting in the **SPECIMEN** category column.
- A.5. Select **FILE>EXPORT TABLE FOR CODIS**. The export dialog box opens:
 - A.5.1 Select **LOOK IN**: File where exported table is
 - A.5.2 Export files as CMF 3.2 (.xml)
 - A.5.3 Source: NV0160000
 - A.5.4 Destination: NV0160000
 - A.5.5 File name: assign file name
- A.6. Click **EXPORT**.
- A.7. Open the **CMF** file in the notebook program.
- A.8. Replace *GMIDX* (or other user name) with a valid CODIS user (i.e. jnaranjo) by using the **FIND/REPLACE** function.

Document ID	Revision	Approval	Date Published
1811	8	Supervising Criminalist - Biology	12/28/2018 11:25:01 AM